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**ONYCHOMYCOSIS TREATMENT - WHITENING AND STRENGTHENING ACTION**  
AGAINST NAIL FUNGUS, YELLOWING, FLAKING AND NAIL LOSS

**INDICATIONS**

Many factors may cause **weakening of the keratin layer** of the nails, increasing the risk of aggression by bacteria and fungi with the resulting unpleasant and unsightly manifestations (onychomycosis, flaking, yellowing and nail loss).

Onychomycosis is encouraged by various situations, such as:

Tight and non-breathing shoes (e.g. sports shoes, safety shoes, work shoes); swimming pools, gyms, changing rooms in sports environments; bumping the nails; reduced peripheral circulation conditions because of vascular problems, diabetes, psoriasis or other.

Onychomycosis may strike one or more nails on both hands and feet, even if it more frequently strikes toenails.

DERMOVITAMINA MICOBLOCK® 3 in 1 ONYCHOMYCOSIS NAIL SOLUTION is not a medicine. It is a medical device for treating and preventing nail mycosis, or in any case a nail suffering from onychomycosis phenomena.

Treats and prevents the nail from yellowing, flaking and falling off. DERMOVITAMINA MICOBLOCK® 3 in 1 ONYCHOMYCOSIS NAIL SOLUTION forms a superficial protective layer that mechanically prevents exchange with the outside environment and deposits acidifying substances on the nail in order to create a hostile environment for the proliferation of fungi that can attack the nail bed. The presence of urea softens the surface of the nail, improves adhesion of the film and penetration of the acidifying substances into the upper layers of the nail.

It is recommended that you speak to your pharmacist or doctor to learn the exact reason for the phenomenon. It is also advisable to consult a doctor if the bottom part of the nail appears unhealthy, the nail is black without having suffered trauma or you suffer from diabetes or immunodeficiency.

**The formula is enriched with substances that strengthen the nail making it less fragile; this also makes the product beneficial where nails are affected by psoriasis or damaged by bumps and traumas and therefore more subject to mycotic contamination. In addition, the formula allows improving nail aesthetics right from the first applications, combating yellowing, flaking and nail loss.**

**\* METAL TESTED**

Nickel < 0,2 ppm  
Chrome < 0,2 ppm  
Cobalt < 0,1 ppm

**INSTRUCTIONS FOR USE**

DERMOVITAMINA MICOBLOCK® 3 in 1 ONYCHOMYCOSIS NAIL SOLUTION can be applied very easily every day on clean and dry nails.

**The nails do not need to be filed before applying DERMOVITAMINA MICOBLOCK® 3 in 1 ONYCHOMYCOSIS NAIL SOLUTION.**

The brush allows optimal application of the product on the whole surface of the nail, on the surrounding skin and on the under part of the free edge of the nail. Invisible and odourless, the product dries very quickly, integrating rapidly with the nail structure.

For prevention, it is advisable to apply the product on the nail affected by mycosis as well as on the nearby nails.

For the first month: apply the product twice a day.

For the subsequent months and until complete regrowth of the healthy nail: apply the product once a day avoiding to wet the nail involved for at least one hour after application.

DERMOVITAMINA MICOBLOCK® 3 in 1 ONYCHOMYCOSIS NAIL SOLUTION is not a nail varnish. If you wish, you can apply a cosmetic nail varnish on its protective film (at least one hour after application). Before reapplying remove any cosmetic nail varnish to allow the product to penetrate.

**TREATMENT DURATION**

The results of the treatment with DERMOVITAMINA MICOBLOCK® 3 in 1 ONYCHOMYCOSIS NAIL SOLUTION can be appreciated with the regrowth of the nail. The treatment of the affected nails has to be continued for several consecutive months to allow the physiological regrowth of the healthy nail, eliminating the damaged part. If symptoms persist and there are no visible improvements in regrowth of the nail, consult your doctor.

**PRACTICAL TIPS**

- wash your hands and feet and thoroughly dry them
- apply the product on well dried nails
- do not share your towels with others
- preferably wear cotton socks and leather shoes
- avoid too tight shoes
- wear slippers in changing rooms, swimming pools, saunas, showers and communal areas as a preventive measure and to avoid spreading bacterial infection

**WARNINGS**

For topical application on the nail - Do not apply on lesions - Do not swallow - Close well after use - Do not use in case of hype sensitivity or allergy to one or more of the product components - Keep out of the reach of children - Do not use after the expiry date indicated on the package - It is recommended not to use the product on children under three years of age - Easily flammable - Irritant to the eyes - Do not use if the package is open or damaged - If used simultaneously with other pharmacological therapies, consult your doctor on how to use the product.

**STORAGE**

Store in a cool and dry place, away from flames or heat sources.

**CONTRAINDICATIONS AND SIDE EFFECTS**

There are no known contraindications and side effects. It is recommended not to use the product if you are sensitive to any of its components. In the event of redness or irritation of the skin, stop application and consult your doctor.

**COMPOSITION**

ALCOHOL DENAT., AQUA, CLIMBAZOLE, 3-O-ETHYL ASCORBIC ACID, HYDROXYPROPYLCELLULOSE, LACTIC ACID, PANTHENOL, PENTYLENE GLYCOL, PIROCTONE OLAMINE, PROPYLENE GLYCOL, SODIUM BENZOTRIAZOLYL BUTYLPHENOL SULFONATE, 4-TERPINEOL, UREA

**PACK**

7ml bottle 

With applicator brush and fitted with safety cap

**MEDICAL DEVICE**

Rev. 02.0 - March 2021

**Leaflet:** PAP 22 Paper - Paper collection.

For disposal check the conditions of your municipality

 **Pasquali S.r.l.**  
V.le Belfiore, 41  
50144 Firenze  
ITALY

**LEGENDA**

Do not dispose of in the environment

Do not use if the packaging is not intact

Do not expose to the sun

Read the instructions for use

Warning: Carefully read the warnings

Recommended temperature between 10°C and 30°C